

## 11.4 Measuring process for measuring head RDS- SP25, - VAST XXT TL1 und TL2, RDS VAST XXT TL3 (RDS-CAA)

### 11.4.1 Technical information VAST XXT product family

RDS-SP25,  
RDS-VAST XXT  
- TL1  
- TL2,  
- TL3

	TL1	TL2	TL3
Probe length LT	30-125 mm	125-250 mm	30-150 mm
Maximum laterally protruding probe length	40 mm	40 mm	65 mm
Recommended work area for lateral probes L <sub>c</sub>	50-70 mm	125-175 mm	50-70 mm
Probe diameter	0.3-8 mm	0.3-8 mm	0.3-8 mm
RC-CAA qualification procedure	No	No	Yes
Max. weight of the probe system	10 g.	10 g.	15 g.

